

## Anti-Ampka1/Ampka2 Rabbit Pab For Wb, Ihc-P, Elisa - Q13131 / P54646

Item Number: CM0014971



### Introduction

Rabbit polyclonal antibody targeting AMPK catalytic subunits alpha-1 and alpha-2, validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for metabolic signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | AMPKa1/AMPKa2 Rabbit pAb                                                                                                  |
| Remarks/Alias    | AMPKa1/AMPKa2                                                                                                             |
| Species          | Human                                                                                                                     |
| GeneID (Human)   | 5562/5563                                                                                                                 |
| GeneID           | 5562/5563                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 13-300 of human PRKAA1/PRKAA2 (NP_006242.5) |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, IHC-P, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                         |
| SWISS            | Q13131/P54646                                                                                                             |
| Protein Weight   | 64kDa/65kDa/62kDa                                                                                                         |
| Shipping         | Ice bag                                                                                                                   |