

Anti-Ampka1/Ampka2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13131 / P54646

Item Number: CM0023791



Introduction

Rabbit polyclonal antibody targeting AMPK α 1/AMPK α 2, the catalytic subunits of AMPK. Validated for WB, IF/ICC, ELISA in human, mouse, rat. Detects ~64/65 kDa bands. Ideal for energy metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	AMPK α 1/AMPK α 2 Rabbit pAb
Remarks/Alias	AMPK α 1/AMPK α 2
Species	Human
GeneID (Human)	5562/5563
GeneID	5562/5563
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 5-280 of human AMPK α 1/AMPK α 2 (NP_006242.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13131/P54646
Protein Weight	64kDa/65kDa
Shipping	Ice bag