


Anti-Ampk α 2 Rabbit Polyclonal Antibody For Energy Metabolism Research - P54646

Item Number: CM0015393



Introduction

AMPK α 2 Rabbit pAb validated for WB, IF/ICC, ELISA; reacts with human and rat. Ideal for studies on AMPK signaling, autophagy, and metabolic diseases. - P54646.

[Learn More](#)

Parameter	Value
Product Name	AMPK α 2 Rabbit pAb
Remarks/Alias	AMPK; AMPK2; PRKAA; AMPK α 2; AMPK α 2
Species	Human
GeneID	5563
Immunogen	Recombinant fusion protein containing amino acids 270-343 of human PRKAA2 (NP_006243.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P54646
Protein Weight	62kDa
Shipping	Ice bag