


Anti-Ampk α 2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P54646

Item Number: CM0013880



Introduction

Validated rabbit polyclonal antibody against human AMPK α 2, the catalytic subunit of the AMPK energy sensor kinase. Designed for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for metabolism and cell signaling studies.

[Learn More](#)

Parameter	Value
Product Name	AMPK α 2 Rabbit pAb
Remarks/Alias	AMPK; AMPK2; PRKAA; AMPK α 2; AMPK α 2
Species	Human
GeneID (Human)	5563
GeneID	5563
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-343 of human AMPK α 2 (NP_006243.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54646
Protein Weight	62kDa
Shipping	Ice bag