

Anti-Ampky1 Polyclonal Antibody For Wb, Elisa - P54619

Item Number: CM0020262



Introduction

Rabbit polyclonal antibody against AMPK γ 1 (PRKAG1), validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying AMPK signaling and energy metabolism. Protein weight: 38kDa.

[Learn More](#)

Parameter	Value
Product Name	AMPK γ 1 Rabbit pAb
Remarks/Alias	AMPKG; AMPK γ 1
Species	Human
GeneID	5571
Immunogen	Synthetic peptide, a synthetic peptide corresponding to a sequence within amino acids 100-200 of human AMPK γ 1 (NP_002724.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54619
Protein Weight	38kDa
Shipping	Ice bag