

Anti-Ampkβ1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y478

Item Number: CM0025442



Introduction

Rabbit polyclonal antibody to AMPKβ1 for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat; detects endogenous AMPKβ1 at approximately 30kDa. Suitable for energy metabolism research.

[Learn More](#)

Parameter	Value
Product Name	AMPKβ1 Rabbit pAb
Remarks/Alias	AMPK; HAMPKb; AMPKβ1
Species	Human
Gene ID (Human)	5564
Gene ID	5564
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human AMPKβ1 (NP_006244.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y478
Protein Weight	30kDa
Shipping	Ice bag