


Anti-Ampkβ1 Polyclonal Antibody (Ko Validated) For Wb, If/Icc, And Elisa - Q9Y478

Item Number: CM0013579



Introduction

Rabbit polyclonal antibody targeting the AMPKβ1 subunit (PRKAB1), validated for WB, IF/ICC, and ELISA applications. Detects human, mouse, and rat AMPKβ1. Suitable for energy metabolism, AMPK signaling, and metabolic disease research.

[Learn More](#)

Parameter	Details
Product Name	[KO Validated] AMPKβ1 Rabbit pAb
Remarks/Alias	AMPK; HAMPKb; β1
Species	Human
GeneID (Human)	5564
GeneID	5564
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human AMPKβ1 (NP_006244.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y478
Protein Weight	30kDa
Shipping	Ice bag