


Anti-Ampkβ1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y478

Item Number: CM0012341



Introduction

Recombinant rabbit monoclonal antibody targeting AMPKβ1 (PRKAB1) for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat AMPKβ1. Ideal for energy metabolism and AMPK signaling research.

[Learn More](#)

Parameter	Value
Product Name	AMPKβ1 Rabbit mAb
Remarks/Alias	AMPK; HAMPKb; AMPKβ1
Species	Human
GeneID (Human)	5564
GeneID	5564
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y478
Protein Weight	30kDa
Shipping	Ice bag