

Anti-Amt Polyclonal Antibody For Wb, If/Icc, Elisa - P48728

Item Number: CM0030822



Introduction

Rabbit polyclonal antibody against human AMT, validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat AMT. Ideal for glycine cleavage system and mitochondrial function studies.

[Learn More](#)

Parameter	Value
Product Name	AMT Rabbit pAb
Remarks/Alias	GCE; NKH; GCST; GCVT; AMT
Species	Human
Gene ID (Human)	275
Gene ID	275
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 107-386 of human AMT (NP_001158184.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48728
Protein Weight	44kDa
Shipping	Ice bag