

Anti-Amy1A Polyclonal Antibody For Wb, Elisa - P04745

Item Number: CM0014071

Introduction

AMY1A rabbit polyclonal antibody validated for WB and ELISA with cross-reactivity to mouse. Designed for human alpha-amylase detection in diagnostic and research settings.

[Learn More](#)

Parameter	Value
Product Name	AMY1A Rabbit pAb
Remarks/Alias	AMY1; AMY1A
Species	Human
GeneID (Human)	276
GeneID	276
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-200 of human AMY1A (NP_004029.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P04745
Protein Weight	57kDa
Shipping	Ice bag