


Anti-Amy2A Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P04746

Item Number: CM0025677



Introduction

Rabbit polyclonal antibody targeting human AMY2A (pancreatic alpha-amylase). Validated for WB, IF-P, ELISA; cross-reacts with mouse and rat. Ideal for carbohydrate metabolism and digestive enzyme studies.

[Learn More](#)

Parameter	Value
Product Name	AMY2A Rabbit pAb
Remarks/Alias	PA; AMY2; AMY2A
Species	Human
GeneID (Human)	279
GeneID	279
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-511 of human AMY2A (NP_000690.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04746
Protein Weight	58kDa
Shipping	Ice bag