

Anti-Amylin/Dap Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P10997

Item Number: CM0007301



Introduction

Rabbit monoclonal antibody targeting human Amylin/DAP (IAPP). Suitable for WB, IHC-P, IF/ICC, and ELISA applications. Cross-reacts with mouse and rat. Ideal for metabolic disease and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	Amylin/DAP Rabbit mAb
Remarks/Alias	DAP; IAP; Amylin/DAP
Species	Human
Gene ID	3375 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-70 of human Amylin/DAP (NP_000406.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10997
Protein Weight	10kDa
Shipping	Ice bag