


Anti-Anapc11 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nyg5

Item Number: CM0013242



Introduction

Rabbit polyclonal ANAPC11 antibody validated in WB, IHC-P, IF/ICC, and ELISA. Detects human and mouse ANAPC11 (APC11, 10 kDa, Q9NYG5). Ideal for cell cycle and ubiquitination studies.

[Learn More](#)

Parameter	Value
Product Name	ANAPC11 Rabbit pAb
Remarks/Alias	APC11; Apc11p; HSPC214; ANAPC11
Species	Human
GeneID (Human)	51529
GeneID	51529
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-84 of human ANAPC11 (NP_001002245.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NYG5
Protein Weight	10kDa
Shipping	Ice bag