


## Anti-Anapc7 Polyclonal Antibody For Wb, Elisa - Q9Ujx3

Item Number: CM0022931



### Introduction

Rabbit polyclonal antibody against human ANAPC7 (APC7, APC/C subunit). Suitable for WB, ELISA; cross-reacts with human and mouse. Use in cell cycle and ubiquitination studies.

[Learn More](#)

| Parameter        | Value                                                                                                   |
|------------------|---------------------------------------------------------------------------------------------------------|
| Product Name     | ANAPC7 Rabbit pAb                                                                                       |
| Remarks/Alias    | APC7; FERBON; ANAPC7                                                                                    |
| Species          | Human                                                                                                   |
| GeneID           | 51434                                                                                                   |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 50-150 of human ANAPC7 (NP_057322.2) |
| Source           | Rabbit                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                   |
| Application      | WB, ELISA                                                                                               |
| Cross Reactivity | Human, Mouse                                                                                            |
| SWISS            | Q9UJX3                                                                                                  |
| Protein Weight   | 63kDa                                                                                                   |
| Shipping         | Ice bag                                                                                                 |