


Anti-Androgen Receptor Monoclonal Antibody For Wb, Ip, Ihc-P, Mihc, If-P, Elisa - P10275

Item Number: CM0002775



Introduction

Rabbit monoclonal antibody targeting human Androgen Receptor (P10275), reactive with human, mouse, and rat. Suitable for WB, IP, IHC-P, IF-P, mIHC, and ELISA applications in androgen signaling and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	Androgen Receptor Rabbit mAb
Remarks/Alias	KD; AIS; AR8; TFM; DHTR; SBMA; HYSPI; NR3C4; SMAX1; HUMARA; Androgen Receptor
Species	Human
GeneID (Human)	367
GeneID	367
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, mIHC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10275
Protein Weight	99kDa
Shipping	Ice bag