


Anti-Androgen Receptor Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P10275

Item Number: CM0023684



Introduction

Rabbit polyclonal anti-Androgen Receptor (AR) antibody for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Studies androgen signaling and nuclear receptor function. SWISS: P10275.

[Learn More](#)

Parameter	Value
Product Name	Androgen Receptor Rabbit pAb
Remarks/Alias	KD; AIS; AR8; TFM; DHTR; SBMA; HYSPI; NR3C4; SMAX1; HUMARA
Species	Human
GeneID (Human)	367
GeneID	367
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Androgen Receptor (NP_000035.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10275
Protein Weight	99kDa
Shipping	Ice bag