


Anti-Ang (Angiogenin) Rabbit Polyclonal Antibody For Wb And Elisa - P03950

Item Number: CM0014199



Introduction

Rabbit polyclonal antibody against human Angiogenin (ANG, RNase5). Suitable for WB and ELISA, cross-reacts with rat. Ideal for studies on angiogenesis, stress response, and rRNA transcription.

[Learn More](#)

Parameter	Value
Product Name	ANG Rabbit pAb
Remarks/Alias	ALS9; RAA1; HEL168; RNASE4; RNASE5; ANG
Species	Human
GeneID (Human)	283
GeneID	283
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-147 of human ANG (NP_001136.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P03950
Protein Weight	17kDa
Shipping	Ice bag