


Anti-Ang Polyclonal Antibody For Wb, If/Icc, Elisa - P03950

Item Number: CM0015458



Introduction

Rabbit polyclonal antibody against human angiogenin (ANG/RNASE5), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for angiogenesis, stress response, and cell proliferation studies.

[Learn More](#)

Parameter	Value
Product Name	ANG Rabbit pAb
Remarks/Alias	ALS9; RAA1; HEL168; RNASE4; RNASE5; ANG
Species	Human
GeneID (Human)	283
GeneID	283
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-147 of human ANG (NP_001136.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P03950
Protein Weight	17kDa
Shipping	Ice bag