


Anti-Angiotensinogen Rabbit Monoclonal Antibody For Wb, Elisa - P01019

Item Number: CM0002299



Introduction

Rabbit monoclonal antibody targeting human angiotensinogen (AGT, Serpin A8). Validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for renin-angiotensin system research in hypertension and cardiovascular studies.

[Learn More](#)

Parameter	Value
Product Name	Angiotensinogen Rabbit mAb
Remarks/Alias	ANHU; hFLT1; SERPINA8; Angiotensinogen
Species	Human
Gene ID (Human)	183
Gene ID	183
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Angiotensinogen (P01019)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01019
Protein Weight	53kDa
Shipping	Ice bag