


Anti-Angiotensinogen Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P01019

Item Number: CM0016210



Introduction

Rabbit polyclonal antibody against human Angiotensinogen (AGT/SERPINA8), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for RAS studies.

[Learn More](#)

Parameter	Specification
Product Name	Angiotensinogen Rabbit pAb
Remarks/Alias	ANHU; hFLT1; SERPINA8; Angiotensinogen
Species	Human
GeneID (Human)	183
GeneID	183
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-485 of human Angiotensinogen (NP_000020.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01019
Protein Weight	53kDa
Shipping	Ice bag