


Anti-Angpt1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q15389

Item Number: CM0012106



Introduction

Anti-ANGPT1 rabbit monoclonal antibody validated for Western blot, IHC-P, and ELISA. Recognizes human, mouse, rat ANGPT1. Ideal for angiogenesis, vascular development, and endothelial biology research.

[Learn More](#)

Parameter	Value
Product Name	ANGPT1 Rabbit mAb
Remarks/Alias	AGP1; AGPT; ANG1; HAE5; AGPT-1; ANGPT1
Species	Human
Gene ID (Human)	284
Gene ID	284
Immunogen	Synthetic peptide (aa 200-300 of human ANGPT1, Q15389)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15389
Protein Weight	58k Da
Shipping	Ice bag