


Anti-Angpt1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15389

Item Number: CM0028210



Introduction

Rabbit polyclonal antibody against Angiopoietin-1 (ANGPT1), validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for angiogenesis and vascular stability studies.

[Learn More](#)

Parameter	Value
Product Name	ANGPT1 Rabbit pAb
Remarks/Alias	AGP1; AGPT; ANG1; HAE5; AGPT-1; ANGPT1
Species	Human
GeneID (Human)	284
GeneID	284
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 277-498 of human ANGPT1 (NP_001137.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15389
Protein Weight	58kDa
Shipping	Ice bag