

Anti-Angpt2 Polyclonal Antibody For Wb, If/Icc, Elisa - 015123

Item Number: CM0013771



Introduction

Rabbit polyclonal antibody against human ANGPT2, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Targets angiogenesis regulator angiopoietin-2 for vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	ANGPT2 Rabbit pAb
Remarks/Alias	ANG2, AGPT2, LMPHM10, ANGPT2
Species	Human
GeneID (Human)	285
GeneID	285
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 197-496 of human ANGPT2 (NP_001138.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	015123
Protein Weight	57kDa
Shipping	Ice bag