


Anti-Angpt2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O15123

Item Number: CM0003023



Introduction

Rabbit monoclonal antibody against human ANGPT2 (angiopoietin-2). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Recognizes a 57kDa protein. Useful for angiogenesis and lymphangiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	ANGPT2 Rabbit mAb
Remarks/Alias	ANG2; AGPT2; LMPHM10; ANGPT2
Species	Human
Gene ID (Human)	285
Gene ID	285
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human ANGPT2 (O15123).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15123
Protein Weight	57kDa
Shipping	Ice bag