


Anti-Angptl3 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Y5C1

Item Number: CM0009856



Introduction

Anti-ANGPTL3 rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat ANGPTL3. Suitable for lipid metabolism, angiogenesis, and adipocyte biology research.

[Learn More](#)

Parameter	Value
Product Name	ANGPTL3 Rabbit mAb
Remarks/Alias	ANL3; ANG-5; FHBL2; ANGPT5; ANGPTL3
Species	Human
Gene ID (Human)	27329
Gene ID	27329
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ANGPTL3 (Q9Y5C1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5C1
Protein Weight	54kDa
Shipping	Ice bag