

Anti-Angptl5 Rabbit Polyclonal Antibody For Wb, Elisa - Q86Xs5

Item Number: CM0018540

Introduction

ANGPTL5 Rabbit pAb validated for WB and ELISA. Reacts with human ANGPTL5, a secreted glycoprotein mainly expressed in adult heart. Ideal for cardiovascular research and biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	ANGPTL5 Rabbit pAb
Remarks/Alias	ANGPTL5
Species	Human
GeneID (human)	253935
GeneID	253935
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-388 of human ANGPTL5 (NP_835228.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86XS5
Protein Weight	44kDa
Shipping	Ice bag