

Anti-Angptl8 Polyclonal Antibody For Wb, Elisa - Q6Uxh0

Item Number: CM0021527

Introduction

Rabbit polyclonal antibody against human ANGPTL8 (betatrophin), validated for WB and ELISA. Detects human, mouse, and rat ANGPTL8 (~22 kDa). Key tool for studying lipid metabolism, triglyceride regulation, and diabetes.

[Learn More](#)

Parameter	Specification
Product Name	ANGPTL8 Rabbit pAb
Remarks/Alias	RIFL; TD26; PRO1185; PVPA599; C19orf80; ANGPTL8
Species	Human
GeneID (Human)	55908
GeneID	55908
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-198 of human ANGPTL8 (NP_061157.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6UXH0
Protein Weight	22kDa
Shipping	Ice bag