


Anti-Ank3 Polyclonal Antibody For If/Icc, Elisa - Q12955

Item Number: CM0021412



Introduction

Rabbit polyclonal antibody targeting human ANK3/Ankyrin-G, validated for IF/ICC and ELISA, cross-reacts with human and mouse. Ideal for neuroscience, muscle, and ion channel research.

[Learn More](#)

Parameter	Value
Product Name	ANK3 Rabbit pAb
Remarks/Alias	MRT37; ANKYRIN-G; ANK3
Species	Human
GeneID (Human)	288
GeneID	288
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q12955
Protein Weight	480kDa
Shipping	Ice bag