

Anti-Ankh Rabbit Polyclonal Antibody For Wb, Elisa - Q9Hcj1

Item Number: CM0028852

Introduction

High-quality Anti-ANKH rabbit polyclonal antibody validated for WB and ELISA. Targets human mineralization regulator ANKH (Q9HCJ1), a key transporter in bone and joint calcification studies.

[Learn More](#)

Parameter	Value
Product Name	ANKH Rabbit pAb
Remarks/Alias	ANK; CMDJ; HANK; MANK; CCAL2; CPPDD; SLC62A1; ANKH
Species	Human
Gene ID (Human)	56172
Gene ID	56172
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-320 of human ANKH (NP_473368.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HCJ1
Protein Weight	54kDa
Shipping	Ice bag