

## Anti-Ankhd1 Polyclonal Antibody For Wb, Elisa - Q8lwz3

Item Number: CM0023143

### Introduction

Anti-ANKHD1 rabbit polyclonal antibody for detection of human ANKHD1 (MASK) in Western blot and ELISA. This scaffolding protein is associated with leukemia cell phenotype and apoptosis regulation. Ideal for research.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | ANKHD1 Rabbit pAb                                                                                                     |
| Remarks/Alias    | MASK; MASK1; VBARP; PP2500; ANKHD1                                                                                    |
| Species          | Human                                                                                                                 |
| GeneID (human)   | 54882                                                                                                                 |
| GeneID           | 54882                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1620-1800 of human ANKHD1 (NP_060217.1) |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, ELISA                                                                                                             |
| Cross Reactivity | Human                                                                                                                 |
| SWISS            | Q8IWX3                                                                                                                |
| Protein Weight   | 269kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                               |