

Anti-Ankhd1 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Iwz3

Item Number: CM0027684

Introduction

Rabbit polyclonal antibody against human ANKHD1, validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for leukemia cell biology, apoptosis, and scaffolding protein research.

[Learn More](#)

Parameter	Value
Product Name	ANKHD1 Rabbit pAb
Remarks/Alias	MASK; MASK1; VBARP; PP2500; ANKHD1
Species	Human
GenelD (human)	54882
GenelD	54882
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8IWZ3
Protein Weight	269kDa
Shipping	Ice bag