

Anti-Ankle2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q86XL3

Item Number: CM0021378

Introduction

High-quality ANKLE2 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA; detects human, mouse, rat. Essential tool for studying mitotic nuclear envelope reassembly and brain development.

[Learn More](#)

Parameter	Value
Product Name	ANKLE2 Rabbit pAb
Remarks/Alias	Lem4; LEMD7; MCPH16; KIAA0692; ANKLE2
Species	Human
Gene ID (Human)	23141
Gene ID	23141
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86XL3
Protein Weight	104kDa
Shipping	Ice bag