

Anti-Ankrd1 Monoclonal Antibody For Wb, Elisa - Q15327

Item Number: CM0005691

Introduction

Rabbit anti-ANKRD1 monoclonal antibody validated for WB and ELISA. Reacts with human and mouse ANKRD1 (36 kDa), a cardiac ankyrin repeat protein involved in endothelial activation and cardiac gene regulation. Ideal for IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	ANKRD1 Rabbit mAb
Remarks/Alias	ALRP; CARP; C-193; CVARP; MCARP; bA320F15.2; ANKRD1
Species	Human
GeneID (Human)	27063
GeneID	27063
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human ANKRD1 (NP_055206.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15327
Protein Weight	36kDa
Shipping	Ice bag