

Anti-Ankrd1 Rabbit Polyclonal Antibody For Wb And Elisa - Q15327

Item Number: CM0026991

Introduction

Rabbit polyclonal antibody targeting human ANKRD1 (cardiac ankyrin repeat protein). Suitable for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for cardiovascular and endothelial research.

[Learn More](#)

Parameter	Value
Product Name	ANKRD1 Rabbit pAb
Remarks/Alias	ALRP; CARP; C-193; CVARP; MCARP; bA320F15.2; ANKRD1
Species	Human
Gene ID (Human)	27063
Gene ID	27063
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human ANKRD1 (NP_055206.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15327
Protein Weight	36 kDa
Shipping	Ice bag