

Anti-Ankrd1 Rabbit Polymab® For Wb And Elisa - Q15327

Item Number: CM0006709

Introduction

High-quality Anti-ANKRD1 Rabbit PolymAb® validated for WB and ELISA. Detects human, mouse, rat ANKRD1 (CARP), a cardiac transcription factor and endothelial activation marker. Suitable for apoptosis and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	ANKRD1 Rabbit PolymAb®
Remarks/Alias	ALRP; CARP; C-193; CVARP; MCARP; bA320F15.2
Species	Human
GeneID (Human)	27063
GeneID	27063
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human ANKRD1 (NP_055206.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15327
Protein Weight	36kDa
Shipping	Ice bag