

Anti-Ankrd28 Rabbit Polyclonal Antibody For Wb, Elisa - O15084

Item Number: CM0028037

Introduction

High-quality rabbit polyclonal antibody against human ANKRD28 (PITK/FAP79), validated for WB and ELISA. Cross-reacts with mouse. Ideal for studying PP6-mediated dephosphorylation and cell migration signaling.

[Learn More](#)

Parameter	Value
Product Name	ANKRD28 Rabbit pAb
Remarks/Alias	PITK; FAP79; CFAP79; PPP1R65; ANKRD28
Species	Human
GeneID (human)	23243
GeneID	23243
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ANKRD28 (NP_056014.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15084
Protein Weight	113kDa
Shipping	Ice bag