

Anti-Ankrd46 Rabbit Polyclonal Antibody - Q86W74

Item Number: CM0021211

Introduction

High-quality rabbit polyclonal antibody against ANKRD46, suitable for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for membrane protein and ankyrin repeat research.

[Learn More](#)

Parameter	Value
Product Name	ANKRD46 Rabbit pAb
Remarks/Alias	ANK-S; GENX-115279; ANKRD46
Species	Human
GenelD (human)	157567
GenelD	157567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-232 of human ANKRD46 (NP_001257308.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q86W74
Protein Weight	25kDa
Shipping	Ice bag