

Anti-Ankrd52 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Nb46

Item Number: CM0024441

Introduction

Rabbit polyclonal antibody targeting human ANKRD52, a putative regulatory subunit of PP6. Suitable for Western blot and ELISA, cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	ANKRD52 Rabbit pAb
Remarks/Alias	ANKRD33; ANKRD52
Species	Human
Gene ID (Human)	283373
Gene ID	283373
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-660 of human ANKRD52 (NP_775866.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NB46
Protein Weight	115kDa
Shipping	Ice bag