


Anti-Anln Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nqw6

Item Number: CM0026085



Introduction

High-quality rabbit polyclonal antibody against human ANLN, validated for WB and ELISA. Suitable for studying cytokinesis, cell division, and podocyte biology. Targets the C-terminal region (aa 845-1124).

[Learn More](#)

Parameter	Value
Product Name	ANLN Rabbit pAb
Remarks/Alias	scra; FSGS8; Scraps; ANLN
Species	Human
Gene ID (Human)	54443
Gene ID	54443
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 845-1124 of human ANLN (NP_061155.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NQW6
Protein Weight	124kDa
Shipping	Ice bag