


Anti-Annexin A11 Monoclonal Antibody For Wb, Elisa - P50995

Item Number: CM0002326



Introduction

Rabbit monoclonal antibody against human Annexin A11, validated for Western blotting (WB) and ELISA. Suitable for human and mouse species. Detects a band at 54 kDa. Ideal for cell cycle and cytokinesis studies.

[Learn More](#)

Parameter	Value
Product Name	Annexin A11 Rabbit mAb
Remarks/Alias	ALS23; ANX11; CAP50; CAP-50; IBMWMA; Annexin A11
Species	Human
GeneID (Human)	311
GeneID	311
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 276-505 of human Annexin A11 (NP_001148.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P50995
Protein Weight	54kDa
Shipping	Ice bag