


Anti-Annexin A11 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P50995

Item Number: CM0026508



Introduction

Rabbit polyclonal antibody against human Annexin A11, validated for WB, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Recombinant immunogen, 54 kDa protein, for cell division and calcium signaling research.

[Learn More](#)

Parameter	Value
Product Name	Annexin A11 Rabbit pAb
Remarks/Alias	ALS23; ANX11; CAP50; CAP-50; IBMWMA; Annexin A11
Species	Human
Gene ID (Human)	311
Gene ID	311
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 276-505 of human Annexin A11 (NP_001148.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50995
Protein Weight	54kDa
Shipping	Ice bag