


Anti-Annexin A1/Anxa1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P04083

Item Number: CM0004648



Introduction

Rabbit monoclonal antibody against human Annexin A1 (ANXA1), validated for Western Blot, IHC, ICC, IP, and ELISA. Cross-reacts with mouse and rat. Ideal for inflammation and immunity research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Annexin A1/ANXA1 Rabbit mAb
Remarks/Alias	ANX1; LPC1
Species	Human
GeneID (Human)	301
Immunogen	Recombinant fusion protein containing amino acids 2-346 of human Annexin A1/ANXA1 (NP_000691.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04083
Protein Weight	39kDa
Shipping	Ice bag