

Anti-Annexin A2 Polyclonal Antibody For Wb - P07355

Item Number: CM0022126



Introduction

Annexin A2 Rabbit Polyclonal Antibody (p36, ANX2, LIP2) validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Recombinant immunogen. Ideal for studying calcium-regulated membrane binding and endosomal stability.

[Learn More](#)

Parameter	Value
Product Name	Annexin A2 Rabbit pAb
Remarks/Alias	P36; ANX2; LIP2; LPC2; CAL1H; LPC2D; ANX2L4; PAP-IV; HEL-S-270; Annexin A2
Species	Human
Gene ID	302
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-339 of human Annexin A2 (NP_004030.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07355
Protein Weight	39kDa
Shipping	Ice bag