

Anti-Annexin A4/Annexin Iv Rabbit Monoclonal Antibody For Wb, Elisa - P09525

Item Number: CM0011466



Introduction

Rabbit monoclonal antibody against human Annexin A4/Annexin IV (UniProt P09525). Validated for Western Blot (WB) and ELISA. Cross-reacts with human, mouse, rat. Suitable for membrane biology and exocytosis research.

[Learn More](#)

Parameter	Value
Product Name	Annexin A4/Annexin IV Rabbit mAb
Remarks/Alias	ANX4; P32.5; PIG28; PP4-X; ZAP36; PAP-II; HEL-S-274; Annexin A4/Annexin IV
Species	Human
GeneID (Human)	307
GeneID	307
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Annexin A4 (P09525)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09525
Protein Weight	36kDa
Shipping	Ice bag