

Anti-Annexin A4/Annexin Iv Rabbit Polyclonal Antibody For Wb And Elisa - P09525

Item Number: CM0013795



Introduction

Rabbit polyclonal antibody targeting human Annexin A4/Annexin IV, validated for WB and ELISA. Cross-reacts with human and mouse. Ideal for studying calcium-dependent membrane fusion and exocytosis.

[Learn More](#)

Parameter	Value
Product Name	Annexin A4/Annexin IV Rabbit pAb
Remarks/Alias	ANX4; P32.5; PIG28; PP4-X; ZAP36; PAP-II; HEL-S-274; Annexin A4/Annexin IV
Species	Human
Gene ID (Human)	307
Gene ID	307
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 112-321 of human Annexin A4/Annexin IV (NP_001144.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09525
Protein Weight	36kDa
Shipping	Ice bag