

Anti-Annexin A5/Annexin V Polyclonal Antibody For Wb, Ihc-P, Elisa - P08758

Item Number: CM0016098



Introduction

Rabbit polyclonal antibody targeting human Annexin A5 (Annexin V), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat; ~36 kDa. Suited for anticoagulant and blood coagulation studies.

[Learn More](#)

Parameter	Value
Product Name	Annexin A5/Annexin V Rabbit pAb
Remarks/Alias	PP4; ANX5; ENX2; CPB-I; RPRGL3; HEL-S-7; VAC-alph; Annexin A5/Annexin V
Species	Human
Gene ID	308
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human Annexin A5/Annexin V (NP_001145.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08758
Protein Weight	36kDa
Shipping	Ice bag