

Anti-Annexin A5/Annexin V Rabbit Monoclonal Antibody For Wb, Elisa - P08758

Item Number: CM0007856



Introduction

High-quality anti-Annexin A5/Annexin V rabbit monoclonal antibody validated for WB and ELISA in human, mouse, and rat samples. 36kDa molecular weight. Ideal for hemostasis and thrombosis research.

[Learn More](#)

Parameter	Value
Product Name	Annexin A5/Annexin V Rabbit mAb
Remarks/Alias	PP4; ANX5; ENX2; CPB-I; RPRGL3; HEL-S-7; VAC-alph; Annexin A5/Annexin V
Species	Human
GeneID (Human)	308
GeneID	308
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 221-320 of human Annexin A5 (NP_001145.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08758
Protein Weight	36kDa
Shipping	Ice bag