

Anti-Annexin A5/Annexin V Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P08758

Item Number: CM0021507



Introduction

High-quality rabbit polyclonal antibody targeting Annexin A5/Annexin V, validated for WB, IF/ICC, ELISA. Reacts with human and mouse. Ideal for studying blood coagulation and hemostasis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Annexin A5/Annexin V Rabbit pAb
Remarks/Alias	PP4; ANX5; ENX2; CPB-I; RPRGL3; HEL-S-7; VAC-alpha; V
Species	Human
GeneID (Human)	308
GeneID	308
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human Annexin A5/Annexin V (NP_001145.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08758
Protein Weight	36kDa
Shipping	Ice bag