


Anti-Annexin A6/Annexin Vi Monoclonal Antibody For Wb, Ihc-P, Elisa - P08133

Item Number: CM0007367



Introduction

Rabbit monoclonal antibody against Annexin A6/Annexin VI, validated for WB, IHC-P, ELISA in human, mouse, rat. 76kDa protein, immunogen aa 324-673. Ideal for calcium signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Annexin A6/Annexin VI Rabbit mAb
Remarks/Alias	p68; p70; ANX6; CBP68; CPB-II; Annexin A6/Annexin VI
Species	Human
GeneID (Human)	309
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08133
Protein Weight	76kDa
Shipping	Ice bag