


Anti-Annexin A6/Annexin Vi Polyclonal Antibody For Wb, Elisa - P08133

Item Number: CM0024516



Introduction

Anti-Annexin A6/Annexin VI Rabbit Polyclonal Antibody, validated for Western blot (WB) and ELISA. Reacts with human and mouse. Targets calcium-binding protein involved in membrane trafficking and CD21 association. Suitable for cell biology research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Annexin A6/Annexin VI Rabbit pAb
Remarks/Alias	p68; p70; ANX6; CBP68; CPB-II; VI
Species	Human
GeneID (Human)	309
GeneID	309
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Annexin A6/Annexin VI (NP_001146.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08133
Protein Weight	76kDa
Shipping	Ice bag